


# ORTHODONTIC ANCHOR SCREW




[bioinfinityorthodontics.com](http://bioinfinityorthodontics.com)






## BRACKET HEAD SCREW



Code	Diameter	L	T	TC
BOSB1406-A	Ø 1.4 mm	6 mm	4.5 mm	1.5 mm
BOSB1408-A		8 mm	6.5 mm	
BOSB1410-A		10 mm	8.5 mm	
BOSB1606-A	Ø 1.6 mm	6 mm	4.5 mm	
BOSB1608-A		8 mm	6.5 mm	
BOSB1610-A		10 mm	8.5 mm	
BOSB1612-A		12 mm	10.5 mm	
BOSB1806-A	Ø 1.8 mm	6 mm	4.5 mm	
BOSB1808-A		8 mm	6.5 mm	
BOSB1810-A		10 mm	8.5 mm	
BOSB1812-A		12 mm	10.5 mm	
BOSB2006-A	Ø 2 mm	6 mm	4.5 mm	
BOSB2008-A		8 mm	6.5 mm	
BOSB2010-A		10 mm	8.5 mm	
BOSB2012-A		12 mm	10.5 mm	

BOSB1406-B	Ø 1.4 mm	6 mm	3.5 mm	2.5 mm
BOSB1408-B		8 mm	5.5 mm	
BOSB1410-B		10 mm	7.5 mm	
BOSB1606-B	Ø 1.6 mm	6 mm	3.5 mm	
BOSB1608-B		8 mm	5.5 mm	
BOSB1610-B		10 mm	7.5 mm	
BOSB1612-B		12 mm	9.5 mm	
BOSB1806-B	Ø 1.8 mm	6 mm	3.5 mm	
BOSB1808-B		8 mm	5.5 mm	
BOSB1810-B		10 mm	7.5 mm	
BOSB1812-B		12 mm	9.5 mm	
BOSB2006-B	Ø 2 mm	6 mm	3.5 mm	
BOSB2008-B		8 mm	5.5 mm	
BOSB2010-B		10 mm	7.5 mm	
BOSB2012-B		12 mm	9.5 mm	


## BUTTON HEAD SCREW



Code	Diameter	L	T	TC
BOSS1406-A	Ø 1.4 mm	6 mm	4.5 mm	1.5 mm
BOSS1408-A		8 mm	6.5 mm	
BOSS1410-A		10 mm	8.5 mm	
BOSS1606-A	Ø 1.6 mm	6 mm	4.5 mm	
BOSS1608-A		8 mm	6.5 mm	
BOSS1610-A		10 mm	8.5 mm	
BOSS1612-A		12 mm	10.5 mm	
BOSS1806-A	Ø 1.8 mm	6 mm	4.5 mm	
BOSS1808-A		8 mm	6.5 mm	
BOSS1810-A		10 mm	8.5 mm	
BOSS1812-A		12 mm	10.5 mm	
BOSS2006-A	Ø 2 mm	6 mm	4.5 mm	
BOSS2008-A		8 mm	6.5 mm	
BOSS2010-A		10 mm	8.5 mm	
BOSS2012-A		12 mm	10.5 mm	

BOSS1406-B	Ø 1.4 mm	6 mm	3.5 mm	2.5 mm
BOSS1408-B		8 mm	5.5 mm	
BOSS1410-B		10 mm	7.5 mm	
BOSS1606-B	Ø 1.6 mm	6 mm	3.5 mm	
BOSS1608-B		8 mm	5.5 mm	
BOSS1610-B		10 mm	7.5 mm	
BOSS1612-B		12 mm	9.5 mm	
BOSS1806-B	Ø 1.8 mm	6 mm	3.5 mm	
BOSS1808-B		8 mm	5.5 mm	
BOSS1810-B		10 mm	7.5 mm	
BOSS1812-B		12 mm	9.5 mm	
BOSS2006-B	Ø 2 mm	6 mm	3.5 mm	
BOSS2008-B		8 mm	5.5 mm	
BOSS2010-B		10 mm	7.5 mm	
BOSS2012-B		12 mm	9.5 mm	

## PALATAL SCREW



Code	Diameter	L	T	TC
BOSP2010-C	Ø 2 mm	10 mm	7 mm	3 mm
BOSP2012-C		12 mm	9 mm	3 mm
BOSP2014-D		14 mm	9 mm	5 mm
BOSP2512-C	Ø 2.5 mm	12 mm	9 mm	3 mm
BOSP2514-D		14 mm	9 mm	5 mm
BOSP2516-D		16 mm	11 mm	5 mm

## INSTRUMENTS



BOS-SDH / Screw Driver Handle



BOS-DSSDH / Driver Shaft for Screw Driver Handle\*



BOS-DSCA20 / Driver Shaft for Contra Angle 20 mm\*



BOS-DSCA30 / Driver Shaft for Contra Angle 30 mm\*

\* Single driver shaft for all BioInfinity screws.

STERILE R  
gamma irradiated

do not reuse




<http://ifu.avrupaimplants.com/en>

**BioInfinity**  
orthodontics

  
**BioInfinity**<sup>®</sup>  
orthodontics



[bioinfinityorthodontics.com](http://bioinfinityorthodontics.com)

 Manufacturer:

**avrupa**  
implant

**Avrupa Implant**

İkitelli OSB HESKOP  
M5 Blok No:82 34490  
Basaksehir / Istanbul - TURKEY  
Phone : +90 212 670 13 00 (pbx)  
Fax : +90 212 670 13 10  
[info@bioinfinityimplants.com](mailto:info@bioinfinityimplants.com)

Rev. No: BOBENG-R001

  
9001:2015

  
13485:2016

  
1984